

# Medonic M-series M32 Innovation built on total quality

#### FOR TODAY'S HEMATOLOGY LABS

Full CBC from one drop of blood Maintenance-free shear valve Unique walk-away Autoloader

Make way for Medonic!



Packed with innovation. Uncompromising accuracy.

Outstanding quality. Medonic M-series M32 analyzers occupy

pride-of-place in the very best hematology labs. There's room for one on your bench-top.

M-series M32S with Autoloader is ideal for small to medium-sized hospitals needing walk-away analysis functionality.

#### Sampling to full CBC in about one minute

Three Medonic M-series M32 analyzers include MPA micro-pipette adapter sampling. Based on a simple fi nger prick sample taken directly from the patient, this unique feature transforms Medonic into what is probably the fastest blood cell counter in the world today.

Simply make a fi nger prick, draw blood into the special 20  $\mu$ l micro-capillary tube, slide it into the adaptor and insert in the analyzer. In about one minute later, you see the full blood status on the touch-sensitive display. There's no preparation required, no pre-dilution, no vacuum tubes and no needles.

Since making a fi nger prick is largely painless, MPA is the ideal way to take and analyze blood samples from children. For blood banks, it's the perfect tool for making fast, pre-donation blood testing. It also saves the vein for donation.

#### Accurate results start with precise collection

In hematology analysis, no one can afford to compromise on quality components or accurate results. We certainly don't. Every Medonic M32 analyzer comes equipped with a high-precision shear-valve. No other sampling technique collects blood as precisely as a shear valve. It's a critical part of the measurement that you'll see refl ected in the accuracy of every result you report.

Furthermore, not only is the Medonic shear valve extremely precise, it's also completely maintenance-free. As well as accurate results, you enjoy significantly reduced maintenance costs. The three-year warranty we give on this component further emphasizes the trust we put in its quality and reliability.

#### Perfect walk-away analyzer

Much as you'll enjoy working with your Medonic M32 analyzer, there will be times when you need to leave it to analyze a whole batch of samples by itself.

For these occasions, the M32S with Autoloader is the ideal instrument to have on your bench. As the top-of-the-range model in the M32-series, it's the perfect walk-away analyzer and ideal for many small to medium-sized hospitals. Just pre-load up to 2 x 20 samples and let it do the work!

Medonic M-series M32 analyzers deliver a full CBC from just 20 μl of blood. Perfect for children and blood banks!

Maintenance-free Medonic M-series M32 shear valve secures accurate results and lowers maintenance costs.

### Pleasure-to-use graphic interface

The user interface is undoubtedly the most visible and best appreciated feature of the Medonic M32-series.



Simple-to-use and easy-on-the-eye, its large 7-inch touch-sensitive display features a 800 x 480 pixel



WVGA color screen with a landscape view.



Different set-up views, each optimized for different user situations, facilitate routine use. The most common include a sample view showing 'classical' graphs and all parameters. A second sample view features the much-used reference-bar display that makes clear indication if results are outside of their expected range.

#### Connect and communicate like never before

Better connectivity and improved communication protocols now let laboratories exchange both patient results and quality control data with interested parties. Medonic M-series M32 instruments are fully-equipped to meet new ways of communicating.

For example, they now feature a host of connection possibilities for printers, keyboards, barcode readers, USBs, etc. This better connectivity includes HL7 industrial standard protocol, now available together with the previous XML protocol. LIS connection is via a LAN port. Furthermore, this improved connectivity is very easy to use, as the front-panel USB port demonstrates.

#### Intuitive software speeds and simplifi es

Powerful software keeps Medonic M-series M32 instruments in the front line of automated hematology analysis, simplifying work processes and giving users better control of sample results and patient records. ID fi eld-length has been increased to 50 characters and patient memory extended to 50,000 records, for example. The menu system is both simple-to-understand and easy-to-navigate.

#### **PLT** extended count

Medonic M-series M32 provide a PLT extended count feature for more accurate PLT results in the critical low PLT range. In case of patient sample with severe thrombocytopenia detected the analyzer will activate PLT extended count, counting three times as many platelets as in normal count mode, to be able to give a more accurate PLT result. PLT results analyzed with extended count will be marked with an asterisk (\*).

#### **Part of our Total Quality Concept**

Medonic M-series M32 hematology analyzers are part of a Total Quality Concept that comprises instruments, reagents, quality control materials and service – an unbroken chain of know-how in quality hematology measurements.

Modern user-interface. Its clean, uncluttered design promotes effi cient operation and accurate assessment of results.



Medonic M-series M32 analyzer users enjoy improved data communication thanks to a broad selection of connectivity options.



Improved software simplifi es work processes and gives users better control of sample results and patient records.



## Four models available One is just right for you





Medonic M32B – even the basic



Medonic M32M – five-sample mixer is ideal for doctors' offices and small labs.



Medonic M32C – closed-tube







<b>Model characteristics</b>	M32B	M32M	M32C	M32S
Built-in tube mixer		•		•
Micro pipette adapter (MPA)		•	•	•
Maintenance-free shear valve	•	•	•	•
Pre-dilution mode	•	•	•	•
Cap-piercing device			•	•
Autoloader				•
Tochnical specification M	odonio M. corio	No M22		

#### **Technical specifi cation Medonic M-series M32**

Parameters	22			
	WBC, LYM, MID, GRAN, LYM%, MID%, GRAN%, RBC, MCV, HCT, PLT,			
	MPV, HGB, MCH, MCHC, RDW%, RDW, PCT, PDW%, PDW, P-LCR, P-LCC			
	≥ 60 samples/hour			
	≤ 50 seconds			
Sampling system	Closed shear valve			
Display	7 inch WVGA true color (24-bit) touch screen			
Interface ports	1 USB device/4 USB host/1 LAN port			
Memory capacity	50,000 samples			
QC software	Yes, including Levey-Jennings and Xb			
	2			
	≤ 1.6%			
	≤ 0.8%			
	≤ 0.5%			
	≤ 0.6%			
	≤ 2.9%			
Dimensions (HWD)	395 x 295 x 475 mm (M32B/M32M/M32C) 395 x 340 x 475 mm (M32S)			
Weight	≤ 18 kg (M32B/M32M/M32C)			
	≤ 22 kg (M32S)			